# BRILLIANT™NG

<b>EN</b>	Instructions for use	2
<b>ES</b>	Instrucciones para el uso	5
<b>PT</b>	Instruções de utilização	9
<b>CS</b>	Návod k použití	12
<b>SK</b>	Návod na použitie	15
<b>PL</b>	Instrukcja stosowania	18
<b>HU</b>	Használati utasítás	22
<b>RO</b>	Instrucțiuni de utilizare	25
<b>HR</b>	Uputstvo za upotrebu	28
<b>BG</b>	Инструкции за употреба	31
<b>RU</b>	Инструкция по применению	35

**Definition**

BRILLIANT NG is a syringeable, sculptable, radiopaque, highly filled nano hybrid composite for the restoration of all teeth. The optimally adjusted shade system of BRILLIANT NG provides two Vita shades in one syringe.

BRILLIANT NG is cured in a tooth cavity or on a model with all common halogen- or LED lights (e.g. Coltène® COLTOLUX LED).

**Composition**

Methacrylates  
Dental glass  
Amorphous silica

**Technical Data**

Average filler particle size:	0.6 µm
Range of particle size:	0.01–2.5 µm
Filler content by volume :	65 %
Filler content by weight:	80 %

**Indications**

BRILLIANT NG is indicated for:

- direct filling of class I, II, III, IV, and V cavities (details see below)
- repair of composite/ceramic veneers
- blocking out of undercuts
- adhesive cementation of ceramic and composite restorations

BRILLIANT NG Dentin is specifically indicated for:

- direct filling of opaque teeth
- direct class V fillings

BRILLIANT NG Enamel is specifically indicated for:

- direct filling of all teeth
- more translucent material with the layering technique

BRILLIANT NG Bleach is specifically indicated for:

- restoration of whitened teeth
- restoration of deciduous teeth
- characterization of chalk spots
- masking of discolored teeth

BRILLIANT NG Transparent is specifically indicated for:

- translucent surface layers in all classes
- reconstruction of incisal edges
- shape and color corrections to improve esthetics
- more translucent material with the layering technique

**Contra-indications**

If the site cannot be isolated after enamel etching, during application, and curing of BRILLIANT NG. If allergies exist to the components of Brilliant NG. If oral hygiene is poor.

**Side effects**

No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.

**Interactions with other agents**

Agents which contain eugenol and/or oil of cloves may affect the polymerization of BRILLIANT NG. The use of zinc-oxide-eugenol cements in combination with BRILLIANT NG should therefore be avoided. Discolorations may occur when using cationic mouth rinses, plaque indicators or chlorhexidine.

**Application**

Cleansing the tooth

The tooth to be treated and the adjacent teeth are cleansed using a brush and prophy-paste which does not contain fluoride.

Selection of shade

Select shades before isolation. Ideally, the shade is determined in broad daylight using the BRILLIANT NG Shade Guide. Shades are identical to those of the shade guide only 24 hours after curing.

Isolation of the cavity (rubber dam)

Adequate isolation is absolutely necessary for achieving optimum results. The use of HYGENIC® /Roeko Dental Dam is recommended.

Preparation of the cavity

When preparing the cavity, every effort should be made to use a technique that is gentle to the tooth substance (principle of Adhesive Restoration Technique). Prepare the enamel and dentine using 80 µm preparation diamonds and finish using 25 µm finishing diamonds. Bevelled enamel margins are recommended to increase the area of adhesion between tooth and filling material, thereby optimizing marginal seal.

Pulp Protection

Covering the areas close to the pulp with hard-setting calcium hydroxide cements using the spot technique will provide protection of the pulp against bacterial infiltration.

Bonding system

Use ONE COAT BOND® SL according to the corresponding instructions for use.

Matrix and interdental wedges

Apply a thin matrix band for class II, III and IV restorations. Fix the matrix proximally using interdental wedges.

Application

Dose BRILLIANT NG directly onto the mixing pad, take up in small increments with Coltène® Composite Instruments or any other application instrument, and apply to the cavity.

Working time

BRILLIANT NG is photosensitive, and should not be exposed to light too long, in particular to the operating light or sunlight. If possible, protect exposed BRILLIANT NG material by means of light opaque covering (not blue).

Polymerization

BRILLIANT NG can be cured using all common halogen- or LED lights (e.g. COLTOLUX LED).

Recommended exposure times:

**Layer Thickness: 2 mm**

	< 800 mW/cm <sup>2</sup>	>800 mW/cm <sup>2</sup>
BRILLIANT NG Dentin A1/B1	30s	20s
BRILLIANT NG Dentin A2/B2	30s	20s
BRILLIANT NG Dentin A3/D3	30s	20s
BRILLIANT NG Dentin A3.5/B3	40s	20s
BRILLIANT NG Dentin C2/C3	40s	20s
BRILLIANT NG Dentin A4/C4	40s	20s
BRILLIANT NG Dentin Bleach	30s	20s

BRILLIANT NG Enamel A1/B1	30s	20s
BRILLIANT NG Enamel A2/B2	30s	20s
BRILLIANT NG Enamel A3/D3	30s	20s
BRILLIANT NG Enamel A3.5/B3	40s	20s
BRILLIANT NG Enamel C2/C3	40s	20s
BRILLIANT NG Enamel A4/C4	40s	20s
BRILLIANT NG Enamel Bleach	30s	20s
BRILLIANT NG Transparent	20s	20s

BRILLIANT NG must be cured layer by layer. Do not remove the inhibition layer since this would affect adhesion between the layers. In the event of underexposure, there is a risk that the ideal strength values may not be achieved. Overexposure is not possible. Equipment whose light output is unknown should be checked according to the corresponding instructions for use.

#### Finishing procedures

Finishing consists of three stages: coarse finishing, fine finishing and polishing.

Finishing diamonds (40 µm and 15 µm) are the most versatile and the least destructive. They should be used with very light pressure, with a constant wiping motion and a generous water spray. The ideal speed ranges between 5000 and 15 000 rpm. Flexible diamond files (80 µm, 40 µm and 15 µm) and strips are suitable for approximal surfaces. For high gloss polishing results use DIATECH Composite polisher or any comparable instruments. Fluoridize all treated surfaces after finishing.

#### **Emergency measures**

In case of direct contact with the oral mucosa, rinsing with tap water is

sufficient. In the case of contact with the eyes, rinse thoroughly with water (10 min) and consult physician.

#### **Note**

Only supplied to dentists and dental laboratories or upon their instructions. Keep out of the reach of children! Do not mix different shades on the mixing pad, since this procedure will lead to air entrapment and porosities.

#### **Storage**

Store at 4–23 °C / 39–73 °F.

Avoid exposure to direct sunlight or other heat sources.

#### **Marking**

The expiry date and LOT number are marked on the syringes.

#### **Caution**

Federal law restricts this device to sale by or on the order of a dentist.

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Coltène/Whaledent AG