# PRESIDENT® light body / regular body

## Definition

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# PRESIDENT® light body / regular body

## Definition

Coltène® PRESIDENT light body / regular body are silicone-based impression materials for use in dentistry with automatic mixing device.

## Material type

Polyvinylsiloxane, addition-type silicone elastomer.

### light body

– low viscosity

Base: green / Catalyst: light green

– ISO 4823, Type 3, low consistency

### regular body

– medium viscosity

Base: blue / Catalyst: light blue

– ISO 4823, Type 2, medium consistency

## Indications

light body :

– Wash material for the putty-wash impression technique

– Syringe material for the simultaneous mixing technique

– Impression material for relining

regular body :

– Impression material for the single-phase impression technique

– Wash material for the putty-wash impression technique

– Syringe material for the simultaneous mixing technique

– Double arch impression technique

– Impression material for relining

## Clinical time

### light / regular body

## Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

## Side effects / Interactions

Polyvinylsiloxanes have a very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient, however, they should remain in the mouth not more than twice the recommended setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

## Trays

Rigid stock or individual trays. For perfect adhesion, we recommend to apply a thin layer of Coltène® Adhesive to all trays. Double arch impression technique: rigid stock partial flexible trays.

## Dispensing preparation

1. Pull off the cartridge plug and discard it.
2. Eject a small amount of material directly out of the cartridge onto a paper towel/tissue held in

your hand. Extrude enough material until it is evident that base and catalyst are being ejected. This ensures proper mixing.

3. Immediately wipe the cartridge orifice to remove excess material.
4. Select the neutral color Mixing tip. Attach the tip to the cartridge by aligning the prongs inside the tip with the cartridge orifice. Seat the tip onto the cartridge and twist a ¼ turn clock-wise (90°) to lock it in place.
5. Push the Oral tip if needed firmly onto the end of the microMixing tip.
6. Squeeze the trigger with a smooth, steady force to mix and dispense material. Avoid abrupt movements. Release the trigger to stop the flow.

## Mixing tip

Disinfect after use by wiping the Mixing tip with a disinfectant. Store in horizontal position. Always store the cartridge with the microMixing tip attached. This seals the cartridge until it is used again and prevents contamination. Only remove mixing tip directly before the device is to be used again, check that material is ejected evenly and attach new mixing tip.

## Impression

Insert the tray and press it into position (2–3 s). Hold it without pressure until material is set. Double arch impression technique: after the introduction of the tray, central occlusion must be achieved and maintained until the material has completely set. Material residue needs considerably more time to set at room temperature. Check intraorally if the material is set before removing from the mouth. High temperatures accelerate and low temperatures slow down setting.

**Important:** In order to guarantee good cohesion with the wash material, the primary impression must be carefully cleaned and dried before further use.

Skin secretion, latex gloves and surfaces contaminated by latex gloves may interfere with the setting process of polyvinylsiloxanes. Only touch material and the surfaces of which the impression is taken (teeth, preparations, retraction cords) with thoroughly washed and rinsed hands (wash for 15 s with a detergent, rinse with water for another 15 s) or with vinyl gloves. Compounds containing eugenol or hemostatic substances can also impede perfect setting.

If H<sub>2</sub>O<sub>2</sub> hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with water afterwards in order to avoid bubble formation.

## Disinfection

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

## Fabrication of models

The stone models may be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (at least 7 days). The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in clear water afterwards. All industry-

standard dental stone model materials (i.e. Fuji-rock Dental Stone, Hard Rock Dental Stone) can be used.

## Electro plating

Compatible with all commercially available copper plating and silver plating baths.

## Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or petroleum ether will dissolve the Coltène® Adhesive. These solvents should only be used in a well-ventilated area. The trays can be cleaned and disinfected as usual.

## Shelf life and storage

Coltène® PRESIDENT fulfills the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed containers at temperatures 15–23 °C / 59–73 °F and 50 % relative air humidity. Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

## Marking

The expiry date and LOT number (LOT) are shown on the package and cartridge.

## Technical data ISO 4823:2000

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

## PRESIDENT light body / regular body

Mixing time:	0:15 min
Total working time:	1:45 min
Oral setting time:	2:00 min

## Instructions for use date

06-2014

Only supplied to dentists and dental laboratories or upon their instructions.

## Caution:

Federal law restricts this device to sale by or on the order of a dentist.

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**Coltène/Whaledent AG** 🇨🇭