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PRESIDENT putty / putty soft / fast putty soft

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putty / putty soft 0:30 min
fast putty soft 0:30 min

Mundverweildauer ISO:

putty / putty soft 3:00 min
fast putty soft 2:00 min

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PRESIDENT® putty / putty soft / fast putty soft

Definition

Coltène® PRESIDENT putty / putty soft / fast putty soft are A-silicone-based impression materials for use in dentistry.

Material type

Polyvinylsiloxane, addition-type silicone elastomer, putty material.

Composition

putty

Base: yellow
Catalyst: light grey
– ISO 4823, Type 0, putty consistency

putty soft / fast putty soft

Base: brown
Catalyst: light grey
– ISO 4823, Type 0, putty consistency

Indications

putty / putty soft / fast putty soft

- Primary impression in the putty-wash impression technique
- Tray material for the simultaneous mixing technique
- Tray material for the copperband impression technique
- Impressions for study models, orthodontic model, matrices

Clinical time

PRESIDENT putty / putty soft

PRESIDENT fast putty soft

Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

Side effects / Interactions

Polyvinylsiloxanes have a very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient, however, they should remain in the mouth not more than twice the setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

Trays

Rigid stock trays. For perfect adhesion, we recommend to apply a thin layer of Coltène® Adhesive to all trays.

Dosage

Measure out level scoopfulls of Base and Catalyst. **Attention:** Avoid interchanging the color coded measuring-scoops, each belonging to one of the components, since there is a risk of contamination!

Mixing

1. Take the base material and the catalyst with the fingertips and start mixing.
2. Knead the material with the fingertips until the color is uniform.

Attention: Only use the Base and Catalyst from the same package.

Impression

Insert the tray and press it into position (2–3 s). Hold it without pressure until material is set. Material residue needs considerably more time to set at room temperature. You should therefore check intraorally if the material is set before removing from the mouth. Intensive mixing or high temperatures speed up, low temperatures slow down the setting.

Important: Skin secretion, latex gloves and surfaces contaminated by latex gloves may interfere with the setting process of polyvinylsiloxanes. Only touch material and the surfaces of which the impression is taken (teeth, preparations, retraction cords) with thoroughly washed and rinsed hands (wash for 15 s with a detergent, rinse with warm water for another 15 s) or with vinyl gloves. Compounds containing eugenol or hemostatic substances can also impede perfect setting. If H₂O₂ hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water afterwards in order to avoid bubble formation. In order to guarantee good cohesion with the wash material, the primary impression must be carefully cleaned and dried before further use.

Disinfection

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

Manufacture of models

Die stone models may be poured after 30 min at the earliest. The impression remains dimensionally stable for practically unlimited period of time (at least 7 days). The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in lukewarm clear water afterwards. All industry-standard dental stone model materials (i.e. Fujirock Dental Stone, Hard Rock Dental Stone) can be used.

Electro plating

Compatible with commercially available copper plating and silver plating baths.

Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or light petrol will dissolve the Coltène® Adhesive. These solvents should only be used in a well-ventilated area. The trays can be cleaned and disinfected as usual.

Shelf life and storage

Coltène® PRESIDENT putty / putty soft / fast putty soft fulfill the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed containers at temperatures 15–23 °C / 59–73 °F and 50 % relative air humidity. Impressions should be stored at normal room temperature. Avoid exposure to direct sunlight or other heat sources.

Marking

The expiry date and number (LOT) are shown on the package and cartridge.

Technical data ISO 4823:2000

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

PRESIDENT putty / putty soft / fast putty soft

Mixing time ISO:	
putty / putty soft	0:30 min
fast putty soft	0:30 min
Oral setting time ISO:	
putty / putty soft	3:00 min
fast putty soft	2:00 min
Recovery after deformation:	99.0 %
Strain in compression:	
putty / putty soft	1.5 – 2.2 %
fast putty soft	1.6 %
Dimensional change:	- 0.2 %

Instructions for use date

09-2014

Caution

Federal law restricts this device to sale by or on the order of a dentist.

Only supplied to dentists and dental laboratories or upon their instructions.

Coltène/Whaledent AG

Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel +41 71 757 5300
Fax +41 71 757 5301
info.ch@coltene.com

For MSDS see

www.coltene.com

PRESIDENT® putty / putty soft / fast putty soft

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