VOCO

GB **DE** **FR** **ES** **PT** **IT** **GR** **NL** **DK** **FI** **NO** **SE** **CE** 0482

in accordance with DIN EN ISO 9917-1

Instructions for use

Product description:
Ionobond is a glass ionomer cement for lining and core build-up. **Ionobond** permits reliable lining for composite materials and amalgam. It sets rapidly and is radiopaque.

Fields of application:
For cavity lining and core build-up.

Application:
After preparing the cavity, clean and dry enamel and dentine thoroughly. Do not use any aggressive cleaning agents. After removing the smear layer prevent contamination by saliva. If a composite filling is to be placed, bevel the edges of the enamel and create retentive undercuts for dentine margins. Shake powder before use.

Ionobond should be worked at room temperature (15 °C - 23 °C). A cooled glass plate will prolong the working time. Use 1 drop of liquid to 1 level measuring scoop of powder (hold bottle vertically). This corresponds to a mixing ratio of 3,6 g powder to 1 g liquid. Place powder and liquid onto a glass plate or suitable mixing pad and mix the powder into the liquid, a little at a time to begin with and then in larger quantities. For mixing use an abrasion resistant metal or plastic spatula.

Mixing time	working time	setting time
30 s	2 - 2,5 min	5 - 6 min

The mixing time is approx. 30 s. **Ionobond** can be worked for approx. 1,5 - 2 min. Apply **Ionobond** in a thin, uniform layer to all dentine surfaces. **Ionobond** will set in the mouth within 3 - 4 min. Under amalgam fillings, sealing of the cavity walls is achieved by applying a layer of an appropriate product. The bond with composites is enhanced by creating retentions in the glass ionomer cement. In preparation, the surface is roughened with a fine-grain diamond approx. 4 - 6 min after application. In case the definite treatment is not carried out immediately, fill the whole cavity with **Ionobond** and during the second treatment session, remove the upper layer to insert the definite filling material. After that, bevel and etch the enamel and apply the composite filling. This procedure additionally increases the bond strength between Ionobond and the composite material.

Indications, precautions:
Rinse metal instruments with cold water immediately after use and before the cement residues can set. For very deep cavities (remaining dentine layer less than 1.5 mm) pulpal protection with an appropriate product is required. Foreign substances on the surfaces may reduce the adhesion and should therefore be removed completely (e.g. varnishes, insulating agents, eugenol residues etc.) with the exception of pulp protection. In case of insufficient dentine thickness irritation of the pulp may occur without appropriate pulp protection. Insufficient quantities of powder (incorrect mixing ratio) reduce stability and may potentially also lead to irritations of the pulp. **Ionobond** contains Na, Ca, Al-silicate, polyacrylic acid and parabenes. **Ionobond** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Storage:
Close the bottle thoroughly after use and store at room temperature and protected against direct sun light. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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