

3M™ Filtek™ P60 Posterior Restorative

Rapid cure time and excellent performance.

packable



With its speed and handling, it's way out front among posterior restoratives. And its stylish beauty makes it a winner with patients too.

3M™ Filtek™ P60 Posterior Restorative offers:

Superior Packability – Packs well without sticking to instrument. Excellent marginal adaptation and interproximal contact.

Easy Sculpting – With exceptional resistance to slumping.

Fastest Placement Technique – Takes only 20 seconds to cure an increment of 2.5mm.

Low Polymerization Shrinkage – Less stress on the tooth, less sensitivity potential for the patient.

Excellent Wear Resistance and Strength – Means you can be confident of long-lasting posterior restorations.

It's everything you could ask for in a posterior restorative.

KEY FEATURES

- EXCELLENT HANDLING
- FASTEST PLACEMENT
- LOW SHRINKAGE
- EXCELLENT WEAR RESISTANCE AND STRENGTH



Placement Technique Summary

Adhesive	49	77	62	111
Liner			45	
Composite for 5mm depth	40	70	40	120
Seal			75	
Total time	89	147	222	231
	Filetek™P60	SureFil™	ALERT™	Solitaire™

Time in Seconds

Does not include manipulation time.

It's the ultimate time-saver. 3M™ Filtek™ P60 Posterior Restorative cures in a quick 20 seconds for a 2.5mm increment. That's time you can apply toward your next procedure.

3M Reliability

3M™ Filtek™ P60 Posterior Restorative

Ordering Information

Ordering

3410 3M™ Single Bond Dental Adhesive System



No. 3410 Includes:

- 2 - 3ml phosphoric etchant syringes
- 6ml adhesive
- Brush handle, brush tips, etchant syringe tips, disposable mixing wells, technique guides

- 3M™ Filtek™ P60 Posterior Restorative and the popular Single Bond dental adhesive work together as a system to give you strong restorations you can be sure will last a long time.
- Recent IADR studies give Single Bond dental adhesive excellent ratings for ease of use and high bond strength properties.

3M™ Single Bond Dental Adhesive Refills

3411	3M™ Single Bond Adhesive refill, 6ml
7523	3M™ Scotchbond™ Phosphoric Etchant, 2 - 3ml syringes

7510 3M™ Vitrebond™ Light Cure Glass Ionomer Liner/Base



No. 7510 Includes:

- 9g powder
- 5.5 ml liquid
- Spoon, mixing pads

It combines light-cure convenience with true glass ionomer fluoride release and bond strength.

3M™ Vitrebond™ Refills

7512P	9g powder
7512L	5.5ml liquid
1961P	4 large mixing pads

4720 3M™ Filtek™ P60 Posterior Restorative Introductory Kit



No. 4720 Introductory Kit Includes:

- 1 - 4g syringes each of 3M™ Filtek™ P60 Posterior Restorative in shades A3, B2 and C2
- 1 - 3ml vial of 3M™ Single Bond Dental Adhesive
- 1 - 3ml syringe of 3M™ Scotchbond™ Etchant and syringe tips
- Accessories

3M™ Filtek™ P60 Posterior Refills

4720 A3	1 - 4g syringe, A3 shade
4720 B2	1 - 4g syringe, B2 shade
4720 C2	1 - 4g syringe, C2 shade

3M™ Filtek™ Flow Flowable Restorative



It flows under pressure, but holds its shape and stays in place prior to light cure. It sets a new standard for convenience in anterior and posterior restorations.

Shade	Item Number
A2	8300A2
A3	8300A3
A4	8300A4
B1	8300B1
B2	8300B2
C1	8300C1

Tip Refills	Item Number
Syringe Tips	3700T

Syringe refills contain:
2 - 2g syringes and 20 disposable tips.

If you have questions about this or any other 3M Dental Product, call our technical hotline at 1-800-634-2249. Or, visit us at www.3M.com/dental. 3M Dental manufactures and markets more than 1,300 products used by dentists to improve the oral health of people around the world. From 3M™ Single Bond Dental Adhesive to 3M™ Sof-Lex™ Finishing and Polishing products, you can count on 3M for reliable products that deliver consistent results.

Filtek, Scotchbond, Sof-lex and Vitrebond are trademarks of 3M. SureFil is a trademark of Dentsply International. ALERT is a trademark of Jeneric Pentron. Solitaire is a trademark of Heraeus Kulzer, Inc.



Dental Products

3M Center, Building 275-2SE-03
St. Paul, MN 55144-1000
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario, Canada N6A 4T1
1-800-265-1840 ext.6229



Recycled paper
40% pre-consumer
10% post-consumer

Printed in the U.S.A.

© 3M 2000
70-2009-2190-9